ADVANCED WOUND CARE | BRASIL

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
## Referências

1. Arapoglou V. et al.; Analgesic efficacy of an ibuprofen-releasing foam dressing compared with local best practice for painful exuding wounds. *Journal of Wound Care* (2011) Vol 20, No 7, 319-325.
2. Woo K, Sibbald G, Fogh K, Glynn C, Krasner D, Leaper D, Osterbrink J, Price P, Teot L. Assessment and management of persistent (chronic) and total wound pain. *Int Wound J* 2008; 5:205–215.
3. Product Monograph: Coloplast and Biatain® are registered trademarks of Coloplast A/S and related companies. The information in this product monograph is based on the most up-to-date information found on the subject by searching Pubmed, and hand searching relevant literature up to June 2015. The information concerning indication for use, precautions on Biatain® – Ibu is according to obtained CE mark, Class III Medical Device. The literature analysis has been performed by Dr. Caroline Dowsett, PhD, Ole K.Nielsen, PhD and reviewed by Flemming Wilhelmsen, MD (Coloplast A/S). All rights reserved. Coloplast A/S, Høtvedvej 1, DK-3050 Humlebaek, Denmark.
4. Center of Evidence-Based Medicine. "Evidence-based medicine is the conscientious, explicit and judicious use of current best evidence in making decisions about the care of individual patients"<sup>25</sup>.
5. Jorgensen, B.; Friis, G.J.; Gottrup, F. Pain and quality of life for patients with venous leg ulcers: proof of concept of the efficacy of Biatain-Ibu, a new pain reducing wound dressing. *Wound Repair Regen.* 2006, 14,, 233-239.
6. Domenech R.P. et al.; Effect of an ibuprofen-releasing foam dressing on wound pain: a real life RCT. *Journal of Wound Care* (2008) Vol 17, No 8.
7. Fogh K. et al.; Clinically relevant pain relief with an ibuprofen-releasing foam dressing: Results from a randomised, controlled double blind clinical trial in exuding, painful venous ulcers. *Wound Rep Reg* (2012) Vol 20 815-821.
8. Andersen KE, Franken CPM, Gad P, Larsen AM, Larsen JR, van Neer PAFA, Vuerstaek J, Wuite J, Neumann HAM. A randomized, controlled study to compare effectiveness of two foam dressings in the management of lower leg ulcers. *Ostomy/Wound Manage* 2002; 48:31–41.
9. Glasziou, P.; Goddard, O.; Heneghan, C.; Shepperd, S.; Meats, E. Glossary of terms in Evidence-Based Medicine. Centre for Evidence-Based Medicine, Institute of Health Sciences Old Road Campus, Headington, Oxford, OX3 7LF, United Kingdom. <http://www.cebm.net/glossary.asp>. 25-9-2006.
10. Jorgenson B. et al.; Combined use of an ibuprofen-releasing foam dressing and silver dressing on infected leg ulcers. *Journal of Wound Care* (2008) Vol 17, No 5, 210-4.
11. Gottrup at al.; Reducing wound pain in venous leg ulcers with Biatain IBU: A randomised controlled double-blind clinical investigation on the performance and safety. *Wound Rep Reg* (2008) Vol 16, 615-625.
12. Steffansen, B.; Herping, S.P.K. Novel wound models for characterizing the effects of exudate levels on the controlled release of ibuprofen from foam dressings. European Wound Management Association. Prague, Czech Republic, 2006.
13. Karsten Fogh, MD, DMSc; Maibrít B. Andersen, MSc, PhD; Morten Bischoff-Mikkelsen, MD; Renzo Bause, MD, DMSc; Markus Zutt, MD, DMSc; Stephan Schilling, MD; Jean-Luc Schmutz, MD; Jesús Borbujo, MD; Juan A. Jimenez, RN; Hugues Cartier, MD; Bo Jørgensen, MD. Clinically relevant pain relief with an ibuprofen-releasing foam dressing: Results from a randomized, controlled, double-blind clinical trial in exuding, painful venous leg ulcers. *Wound Repair and Regeneration* 20 (2012) 815-821.
14. Holloway S, Ahmajärvi K, Frescos N, Jenkins S, Oropallo A, Slezáková S, Pokorná A. Holistic Management of Wound-Related Pain. *J Wound Management*, 2024;25(1 Sup1). S1-S84. DOI: 10.35279/jowm2024.25.01.sup01.
15. Liberato SMD, Souza AJGD, Costa IKF, Torres GDV, Vitor AF, Lira ALBDC. A enfermagem no manejo da dor em pessoas com úlcera venosa: revisão integrativa Nursing in the management of pain in people with venous ulcer: integrative review. *Rev Pesqui Cuid É Fundam Online*. 2016 Apr 4;8(2):4109–20.
16. Gethin G, Probst S, Stryja J, Christiansen N, Price P. Evidence for person-centred care in chronic wound care: A systematic review and recommendations for practice. *J Wound Care*. 2020 Sep 1;29(Sup9b): S1–22.

17. Thomas DC, Chui PL, Yahya A, Yap JW. Systematic review of patient education for pressure injury: Evidence to guide practice. *Worldviews Evid Based Nurs*. 2022 Aug;19(4):267–74.
18. Kane FM, Brodie EE, Coull A, Coyne L, Howd A, Milne A, et al. The analgesic effect of odour and music upon dressing change. *Br J Nurs*. 2004 Oct 1;13(Sup4):S4–12.
19. Van Hecke et al<sup>147</sup> identified that pain, discomfort and a lack of valid lifestyle advice by healthcare professionals (including conflicting advice, no specific advice or impossibility of instructions) were identified as primary reasons for non-adherence.
20. Weller et al<sup>148</sup> reviewed the benefits and harms of interventions to help people adhere to compression therapy for VLU. One of the outcomes of the review included pain.
21. Zeleníková R, Vyhlídalová D. Applying honey dressings to non-healing wounds in elderly persons receiving home care. *J Tissue Viability*. 2019 Aug;28(3):139–43.
22. Dubhashi SP, Sindwani RD. A Comparative Study of Honey and Phenytoin Dressings for Chronic Wounds. *Indian J Surg*. 2015 Dec;77(Suppl 3):1209–13.
23. Tsihlikidou A, Govina O, Vasilopoulos G, Kavga A, Vastardi M, Kalemiker-akis I. Intervention for symptom management in patients with malignant fungating wounds – a systematic review. *J BUON Off J Balk Union Oncol*. 2019;24(3):1301–8.
24. Ma X, Zheng C, Hu C, Rahman K, Qin L. The genus *Desmodium* (Fabaceae)-traditional uses in Chinese medicine, phytochemistry and pharmacology. *J Ethnopharmacol*. 2011 Nov;138(2):314–32.
25. Ongarora BG. Recent technological advances in the management of chronic wounds: A literature review. *Health Sci Rep*. 2022 May;5(3):e641.
26. Peplow PV, Chung TY, Baxter GD. Application of low-level laser technologies for pain relief and wound healing: overview of scientific bases. *Phys Ther Rev*. 2010 Aug;15(4):253–85.
27. Geraldo MaGela SaloMél<sup>1</sup>; lydia MaSako Ferreira, TCBC-SP1. Impact of non-adherent Ibuprofen foam dressing in the lives of patients with venous ulcers. *Revista do Colégio Brasileiro de Cirurgiões*. 2017;44(2):116-24.
28. Price P, Fogh K, Glynn C, Krasner D, Osterbrink J, Sibbald RG. Managing painful chronic wounds: the wound pain management model. *Int Wound J* 2007;4(Suppl 1):4-14.
29. R Gary Sibbald, Patricia Koota, Marjorie Fierheller, Kevin Woo. A pilot (real-life) randomised clinical evaluation of a pain-relieving foam dressing: (ibuprofen-foam versus local best practice). *Int Wound J* 2007;4(Suppl. 1):16-23.
30. Rossella Laurano, Monica Boffito, Gianluca Ciardelli, Valeria Chiono. Wound dressing products: A translational investigation from the bench to the market. *Engineered Regeneration* 3 (2022) 182–200.

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